GEM Microvascular Anastomotic COUPLER Device

en	Instructions for use.....	6
fr	Instructions d'utilisation.....	16
de	Gebrauchsanweisung.....	26
it	Istruzioni per l'uso.....	36
es	Instrucciones de uso.....	46
nl	Gebruiksaanwijzing.....	56
da	Betjeningsvejledning.....	66
sv	Bruksanvisning.....	76
tr	Kullanım Talimatları.....	86
el	Οδηγίες χρήσης.....	96
ja	取扱説明書.....	108
zh	使用说明书.....	118
no	Bruksanvisning.....	126
pt	Instruções de utilização.....	136
pl	Instrukcja użycia.....	146
ko	사용 설명서.....	158
is	Notkunarlæðbeiningar.....	168
fi	Käyttöohjeet.....	178

Microvascular Anastomotic COUPLER Device
Dispositif microvasculaire anastomotique COUPLER
COUPLER-Gerät für die mikrovaskuläre Anastomose
Dispositivo microvascolare anastomotico COUPLER
Dispositivo microvascular anastomótico COUPLER
Microvasculair anastomotisch COUPLER-hulpmiddel
Mikrovaskulær anastomotisk COUPLER-enhed
Mikrovaskulär anastomotisk COUPLER-enhet
Mikrovasküler Anastomotik COUPLER Cihazı
Συσκευή μικροαγγειακής αναστόμωσης COUPLER
微小血管吻合用COUPLERデバイス
微血管吻合 COUPLER 器械
Mikrovaskulær anastomotisk COUPLER-enhet
Dispositivo anastomótico microvascular COUPLER
Urządzenie COUPLER do zespozeń mikronaczyniowych
미세혈관 문합기 COUPLER 기기
COUPLER-smáæðasamtengibúnaður
Mikrovaskulaarinen anastomoosin COUPLER-laite


Figure 1

Figure 2


Figure 3


Figure 4


Figure 5a


Figure 5b


Figure 6


Figure 7


Figure 8



Figure 9


Figure 10


Figure 11


Figure 12


Figure 13


Figure 14


Figure 15


Figure 16


Figure 17


Figure 18


Figure 19


Figure 20a


Figure 20b


Figure 21a


Figure 21b


Figure 22



Figure 23

SYMBOL GLOSSARY PER US FD&C ACT:

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO-15223-1*		Manufacturer	Manufacturer	5.1.1
ISO-15223-1		Authorized representative in the European Community/ European Union	Authorized representative in the European Community/ European Union	5.1.2
ISO-15223-1		Use-by Date	Use-by Date	5.1.4
ISO-15223-1		Batch code	Lot number	5.1.5
ISO-15223-1		Catalogue number	Catalogue number	5.1.6
ISO-15223-1		Importer	Indicates the entity importing the medical device into the locale	5.1.8
ISO-15223-1		Country of Manufacture	Identify the country of manufacture of product. The date of manufacture will be adjacent to this symbol	5.1.11
ISO-15223-1		Sterilized using irradiation	Sterilize using irradiation	5.2.4
ISO-15223-1		Do not re-sterilize	Do not re-sterilize	5.2.6
ISO-15223-1		Do not use if package is damaged and consult instructions for use	Do not use if the product sterilization barrier or its packaging is compromised.	5.2.8

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO-15223-1		Single sterile barrier system	Indicates one sterile barrier system	5.2.11
ISO-15223-1		Double sterile barrier system	Indicates two sterile barrier systems	5.2.12
ISO-15223-1		Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed	5.3.7
ISO-15223-1*		Do not re-use	Do not re-use	5.4.2
ISO-15223-1		Consult instructions for use	Consult instructions for use	5.4.3
ISO-15223-1		Patient identification	Indicates the identification data of the patient	5.7.3
ISO-15223-1		Patient information website	Website where a patient may obtain additional information on the medical product	5.7.4
ISO-15223-1		Health care centre or doctor	Address of the health care centre or doctor where medical information about the patient may be found	5.7.5
ISO-15223-1		Date	Indicates the date that information was entered, or a medical procedure took place	5.7.6
ISO-15223-1		Medical device	Indicates the item is a medical device	5.7.7
ISO-15223-1		Unique Device Identifier (UDI)	Indicates a carrier that contains Unique Device Identifier information	5.7.10

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ASTM F2503-23		MR Conditional with conditions	Indicates a demonstrated safety in the MR environment within defined conditions	
			See IFU for symbol definition	
			CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.	
			Indicates the authorized representative in Switzerland	

*ISO 15223-1 Medical devices - Symbols to be used with information to be supplied by the manufacturer. - Part 1: General requirements.

Additional symbols and graphics on the product labeling that are not derived from standards:

Symbol	Symbol Description
	Manufacturer part number
	Internal tracking number
	Serial number
	This product and package is not made with natural rubber latex.
	Content
	Size of the GEM COUPLER Device (inner diameter of the COUPLER rings)
	Microvascular Anastomotic COUPLER Device

DESCRIPTION:

The Synovis MCA GEM Microvascular Anastomotic COUPLER Device (the COUPLER Device) and System (the COUPLER System) have been specifically designed for use in the anastomosis of blood vessels.

The COUPLER Device rings are a permanent implant made of high-density polyethylene and surgical grade stainless steel pins with a combined mass of up to 0.0304 grams. A protective cover and jaw assembly protect the rings and allow for easy loading onto the Anastomotic Instrument. The jaw assembly consists of jaws (Polyetherimide), adapters (Polyetherimide), hinge pin (Stainless Steel), Torsion Spring (Stainless Steel) and holder cover (ABS Plastic). Both the protective cover and jaw assembly are disposable. The COUPLER Device is single-use and available in various sizes.

The COUPLER System consists of a reusable Anastomotic Instrument (surgical-grade stainless steel and titanium), reusable COUPLER Forceps (surgical-grade stainless steel), a reusable Double-Ended Vessel Measuring Gauge (surgical-grade stainless steel), and a Sterilization Tray (anodized aluminum). The COUPLER System is available as a complete system or as single items.

CLINICAL BENEFITS AND PERFORMANCE OBJECTIVES:

Reduced anastomotic and flap ischemic time, reduced rate of thrombosis, vessel size discrepancy apposition correction, and patency rates for venous anastomoses comparable to hand suturing.

Performance Characteristics:

The COUPLER Device provides patent end-to-end and end-to-side anastomoses (with intima-to-intima contact) in veins and arteries of the peripheral vascular system.

ADDITIONAL INFORMATION FOR PATIENT AND HEALTHCARE PROFESSIONALS**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)**

For pertinent safety and clinical information about this product, the Summary of Safety and Clinical Performance intended for healthcare professionals (Part I) and patients (Part II) can be found at the following location: <https://ec.europa.eu/tools/eudamed>.

Until SSCPs are available on the EUDAMED website, you can find the SSCP at <https://sscp.baxter.com>.

Basic UDI-DI:47818000000000000001157R

INTENDED PURPOSE:

The COUPLER Device is intended to be used in the anastomosis of veins and arteries.

INTENDED PATIENT GROUP:

The COUPLER Device is intended for patients who need to undergo anastomosis of veins or arteries for peripheral vascular reconstructive procedures.

Safe use of the COUPLER Device for the anastomosis of growing vessels in children or adolescents has not been established. Not intended for fetal use.

INTENDED USERS:

The COUPLER Device is intended to be used by surgeons but may be prepared by a nurse or a scrub technician.

INDICATIONS FOR USE:

The COUPLER Device is indicated for use with veins and arteries having an outside diameter no smaller than 0.8 mm and no larger than 4.3 mm and a wall thickness of 0.5 mm or less, normally encountered in microsurgical procedures which require vascular reconstitution only in the peripheral vascular system.

CONTRAINDICATIONS:

The COUPLER Device is not indicated for use in patients presenting conditions that would normally preclude microvascular repair with suture technique. Examples of such conditions include, but are not limited to:

- Pre-existing or suspected peripheral vascular disease,
- Ongoing irradiation of the area of reconstruction,
- Clinical infection of the area of reconstruction,
- Anticipated infection due to significant contamination of the area of reconstruction,
- Friability of the vascular tissue due to sclerotic conditions,
- Concurrent diabetes mellitus, or
- Concurrent corticosteroid therapy

UNDESIRABLE SIDE-EFFECTS/ADVERSE EFFECTS

Use of the COUPLER Device involves potential risks normally associated with any implanted device as well as: allergic reaction, delay in revascularization, inflammation at anastomotic site, granuloma, thrombosis /emboli, vessel stasis, stricture, pain, necrosis, fever, chills, infection, septic shock, tissue damage, surgical intervention for device removal, foreign body reaction/erosion, and hemorrhage.

ADVERSE EVENTS REPORTING:

If during the use of this device, or as a result of its use, a serious incident has occurred, please report this incident to the manufacturer, and/or its authorized representative, and/or the competent authority of the Member State and/or relevant regulatory bodies in which the user and/or patient is established.

For reporting a device malfunction or adverse event, or for product inquiries, contact: Synovis MCA sales representative, or email SMCA_Quality@baxter.com.

WARNINGS:

- Failure to use the Vessel Measuring Gauge to approximate the vessel size could result in using a COUPLER Device of an inappropriate size. Using a ring too large for the vessel may result in stressing or tearing of the vessel wall and a compromised anastomosis. Using a ring too small for the vessel may unduly constrict the vessel and lead to thrombosis or ring separation.
- Failure to squeeze the COUPLER jaws with a hemostat or similar instrument prior to ejection of the joined rings may result in an inadequate friction fit and possible ring separation.

Inspect the anastomotic site to ensure that the anastomosis has been satisfactorily completed.

- The COUPLER Device is supplied sterile and is single use only. **Do not resterilize or reuse the COUPLER Device.**
 - Resterilization may compromise the structural integrity of the product which may lead to incomplete anastomosis.
 - Device cannot be reused due to possible structural damage during first use which may lead to incomplete anastomosis.
- Do not use the COUPLER Device if the package has been opened or appears to be damaged or compromised as sterility may be compromised. Failure to observe this warning may result in surgical infection.
- Safe use of the COUPLER Device for the anastomosis of tubular structures other than veins and arteries has not been established.
- Security of an anastomosis utilizing COUPLER Devices that have been approximated, reopened, and then reapproximated has not been demonstrated. When reapproximation of the anastomosis is desired, the vessel should be removed from each ring and a new COUPLER Device utilized.
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray **must be sterilized prior to use.**
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray should be thoroughly inspected before use. Instruments that are damaged and/or in need of repair should not be used.

- When performing an end-to-side anastomosis with this COUPLER device, the lumen of the “side” vessel is narrowed slightly. For this reason, the diameter of the “side” vessel should be larger than that of the “end” vessel when completing such a procedure. The opening made to the side vessel should be approximately the same size as the inside diameter of the COUPLER device being used.

INSTRUCTIONS FOR USE:

These Instructions for Use are designed for proper use of this device. They are not intended to serve as a reference to surgical technique, to supersede institutional protocols or professional clinical judgment regarding patient care.

Read the Instructions for Use before use. No special training is required beyond reading the Instructions for Use.

NOTE: The electronic version of the Instructions for Use can be found at the following location: <https://edocs.baxter.com>.

MRI SAFETY INFORMATION COUPLER RINGS:

It is the responsibility of the clinician to inform the patient that he/she is the recipient of permanent implants which contain metal components (surgical-grade stainless steel pins). The COUPLER Device rings have been evaluated with a 1.5 Tesla and a 3.0 Tesla magnetic field and no change in displacement was observed in each of three orthogonal planes.^{1,3}

The stainless steel pins in the COUPLER Devices are nominally nonferromagnetic. However, the US Food and Drug Administration (FDA) has made recommendations for any medical device implanted which have metallic components to include:

- Documentation in the official medical record of the identity of the implant (manufacturer, model number, lot and serial numbers, and identifying marks, if any).
- Documentation of the technique and results of any magnetic testing performed on the implant or that no such testing was done.
- Patient education regarding the particular implant and recommendation for identifying medical alert card, bracelet, or necklace characterizing the implanted device.²

MRI Information	
Device Name	COUPLER Device
Static Magnetic Field Strength (Bo)	1.5 T or 3.0 T
Maximum Spatial Field Gradient	40 T/m (4,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	There are no Transmit Coil Restrictions
RF Receive Coil Type	Any
Operating Mode	Normal Operating Mode (subject to scan region limitations listed below)
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
Scan Regions	This device can be safely scanned in any region at either 1.5 or 3.0T.
MR Image Artifact	The presence of this implant may produce an image artifact of up to 9 mm from the construct.

END-TO-END ANASTOMOSIS:

Using conventional microsurgical technique, mobilize a minimum of 1 cm of each vessel end. Using vascular clamps, clamp off the vessel(s) and irrigate the vessel openings. The COUPLER Device requires a greater amount of free vessel within the clamps than a conventional suture repair.

1. After gentle dilation, estimate the **outer** diameter of each vessel using the Vessel Measuring Gauge. The circular guides on the gauge **should not** be placed inside the vessel lumen (See Figure 1). If there is a size discrepancy between the two vessels, use the measurement of the smaller vessel to choose the appropriate COUPLER Device. The degree of vessel spasm and the elasticity of the vessel should be considered when choosing the COUPLER Device size to be used.
2. Select the appropriate size COUPLER Device. Both vessel ends should be approximately the same size as the inside diameter of the COUPLER Device being selected.

NOTE: Approximate the vessel to an appropriate COUPLER size using a Double-Ended Vessel Measuring Gauge. In an end-to-end anastomosis, the two vessel ends should be approximately the same size as the inside diameter of the COUPLER Device being used. In an end-to-side anastomosis, the opening made to the side vessel should be approximately the same size as the inside diameter of the COUPLER Device being used. If the end vessel size is outside of the specified range, do not use the COUPLER Device to mate the vessels.

3. Remove the lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact. Remove the lid from the inner tray.

4. Turn the Anastomotic Instrument knob fully counter-clockwise, and then insert the COUPLER Device onto the Anastomotic Instrument. **The matching indicator arrows on the COUPLER Device and the Anastomotic Instrument should be pointing toward each other when loading** (See Figures 2 & 3). **Ensure that an audible click is heard for proper loading.**
5. Remove the COUPLER Device from the protective cover by pulling it firmly away from the Anastomotic Instrument (See Figure 4).
6. Visually inspect to see that both rings are seated at the bottom of the U portion of the jaw and the pins are not bent (See Figures 5a & 5b). If pins are bent, do not attempt to straighten. Instead use a new COUPLER Device.
7. Place the Anastomotic Instrument perpendicular to the vessel(s), with the COUPLER jaw assembly near the two vessel ends. Pull one vessel end through one of the COUPLER rings using microsurgical forceps (See Figure 6).
- 8a. **3.0mm COUPLER Device Size or Smaller:** Take a bite of approximately one to two pin diameters of the vessel wall and intimal lining, evert 90 degrees and impale onto one pin. Proceeding in a triangular fashion, impale the vessel firmly upon every other pin, completing three pins (See Figure 7). Complete vessel placement on the ring by impaling the vessel upon the remaining three intermediate pins (See Figure 8). Ensure that both the vessel wall and the intimal layer are fully impaled upon each pin to reduce the risk of thrombosis. Should the vessel wall tear during impalement, remove the vessel, trim the end, and repeat the procedure. For examples of improper impalement of the vessel see Figure 9.

- 8b. **3.5mm COUPLER Device Size or Larger:** Take a bite of approximately one to two pin diameters of the vessel wall and intimal lining, evert 90 degrees and impale onto the pin situated nearest to the open part of the Jaw Assembly (open end of the U portion of the jaw). Impale the opposite side of the vessel opening to the pin directly across from the initial pin. Next, impale the vessel onto the pins located near the sides of the ring, keeping the vessel as evenly spaced as possible between the four pins (See Figure 15). Continue vessel placement on the ring by impaling the vessel onto the two remaining pins near the open end of the Jaw Assembly. Complete by impaling the vessel onto the last two pins near the bottom of the Jaw Assembly (bottom of the U portion of the jaw); this final step prevents the ring from sliding out of the Jaw Assembly prematurely (See Figure 16). Ensure that both the vessel wall and the intimal lining are fully impaled upon each pin to reduce the risk of thrombosis. Should the vessel wall tear during impalement, remove the vessel, trim the end, and repeat the procedure. For examples of improper impalement of the vessel see Figure 17.
9. Repeat Steps 7 and 8 to impale the other vessel end upon the second COUPLER ring.
10. When both vessel ends have been suitably impaled, visually inspect to ensure that both rings are seated at the bottom of the U portion of the jaw and the pins are not bent (See Figures 5a and 5b). Bring the rings together (See Figures 10 & 11) by turning the Anastomotic Instrument knob clockwise. Turn the knob only until the ejector rod has just started to move the now joined rings.

11. Prior to ejecting the joined rings, gently squeeze the end of the apposed jaws with a small hemostat (See Figure 12) to ensure ring approximation and a tight friction fit. Turn the Anastomotic Instrument knob further clockwise to eject the joined rings.
12. Check the anastomosis under the operating microscope before opening the vascular clamps. Remove the clamps and inspect the anastomotic site to ensure that the anastomosis has been satisfactorily completed (patent vessel without leakage).
13. To remove the jaw assembly turn the Anastomotic Instrument knob fully counterclockwise (See Figure 13). Press the release button, located near the arrow on the Anastomotic Instrument, and remove the jaw assembly (See Figure 14).
14. Rinse the Anastomotic Instrument with water after use.

NOTE: The COUPLER instruments should be cleaned and sterilized in accordance with the MAC instruments system IFU.

END-TO-SIDE ANASTOMOSIS:

Using conventional microsurgical technique, mobilize a minimum of 1 cm of the “end” vessel. Clamp-off the vessel and irrigate the vessel lumen. Mobilize a minimum of 2 cm of the “side” vessel and clamp off the vessel.

1. When performing an end-to-side anastomosis with this COUPLER device, the lumen of the “side” vessel is narrowed slightly. For this reason, the diameter of the “side” vessel should be larger than that of the “end” vessel when completing such a procedure.

2. Estimate the outer diameter of the “end” vessel using the Vessel Measuring Gauge. The circular guides on the gauge **should not** be placed inside the vessel lumen (See Figure 1).
3. Select the appropriate size COUPLER Device.
4. Remove the lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact. Remove the lid from the inner tray.
5. Turn the Anastomotic Instrument knob fully counterclockwise, and then insert the COUPLER Device onto the Anastomotic Instrument. **The matching indicator arrows on the COUPLER Device and the Anastomotic Instrument should be pointing toward each other when loading** (See Figures 2 & 3). **Ensure that an audible click is heard for proper loading.**
6. Remove COUPLER Device from the protective cover by pulling it firmly away from the Anastomotic Instrument (See Figure 4).
7. Visually inspect to see that both rings are seated at the bottom of the U portion of the jaw and the pins are not bent (See Figures 5a and 5b). If pins are bent, do not attempt to straighten. Instead use a new COUPLER Device.
8. Place the Anastomotic Instrument handle perpendicular to the direction of the “end” vessel. Place the “end” vessel upon one ring as described in Steps 7 and 8 for End-to-End Anastomosis directions of the appropriate COUPLER Device size.
9. Create a transverse incision in the “side” vessel of a length no greater than the internal diameter of the COUPLER Device selected. Slide the clamps together slightly to remove tension and to open the incision (See Figure 18). Irrigate the vessel lumen through the created opening.
10. Using microsurgical forceps, grasp the vessel wall and intimal lining near one end of the transverse incision and pull them through the remaining ring. Evert the vessel wall and intimal lining 180 degrees and impale the vessel first upon the pins situated nearest the incision end (See Figure 19).
11. Proceed in a similar manner at the opposite end of the incision, impaling the vessel wall and intimal lining upon the pins situated nearest the incision end (See Figure 20a for COUPLER Device sizes 3.0mm and smaller; See Figure 20b for COUPLER Device sizes 3.5mm and larger). Complete the pinning procedure by everting the vessel upon the remaining pins (See Figure 21a for COUPLER Device sizes 3.0mm and smaller; See Figure 21b for COUPLER Device sizes 3.5mm and larger). Ensure that both the vessel wall and the intimal layer are fully impaled upon each pin.
12. Bring the rings together by turning the Anastomotic Instrument knob clockwise, **ONLY** until the ejector rod has just started to move the now joined rings. Hold the Anastomotic Instrument such that the “end” pinned ring is brought to the “side” pinned ring during the approximation (See Figure 22).

13. Prior to ejecting the joined rings, squeeze the end of the apposed jaws with a small hemostat (See Figure 23) to ensure ring approximation and a tight friction fit. Turn the Anastomotic Instrument knob further clockwise to eject the joined rings.
14. Check the anastomosis under the operating microscope before opening the vascular clamps. Remove the clamps and inspect the anastomotic site to ensure that the anastomosis has been satisfactorily completed (patent vessel without leakage).
15. To remove the jaw assembly turn the Anastomotic Instrument knob fully counterclockwise (See Figure 13). Press the release button, located near the arrow on the Anastomotic Instrument, and remove the jaw assembly (See Figure 14).
16. Rinse the Anastomotic Instrument with water after use.

NOTE: The COUPLER instruments should be cleaned and sterilized in accordance with the MAC instruments system IFU.

SERVICE

For Customer or Technical service, contact:

Phone: +1 205.941.0111 or 1.800.510.3318 (U.S. only)

Fax: +1 205.941.1522

Website: synovismicro.com

ACCESSORIES & PARTS:

ITEM	REF
COUPLER, 1.0MM, 1PK	GEM2750/I
COUPLER, 1.0MM, 6PK	GEM2750
COUPLER, 1.5MM, 1PK	GEM2751/I
COUPLER, 1.5MM, 6PK	GEM2751
COUPLER, 2.0MM, 1PK	GEM2752/I
COUPLER, 2.0MM, 6PK	GEM2752
COUPLER, 2.5MM, 1PK	GEM2753/I
COUPLER, 2.5MM, 6PK	GEM2753
COUPLER, 3.0MM, 1PK	GEM2754/I
COUPLER, 3.0MM, 6PK	GEM2754
COUPLER, 3.5MM, 1PK	GEM2755/I
COUPLER, 3.5MM, 6PK	GEM2755
COUPLER, 4.0MM, 1PK	GEM2756/I
COUPLER, 4.0MM, 6PK	GEM2756
MAC SYSTEM 2 C FORCEPS	GEM2741CC
STERILIZATION TRAY, MAC	GEM2745
ANASTOMOTIC INSTRUMENT, TI	GEM2740
VESSEL MEASURING GAUGE	GEM2749
FORCEPS, COUPLER	GEM4183C

STORAGE CONDITIONS:

The product can be stored at 15°C/59°F to 30°C/86°F.

DISPOSAL METHOD:

Any open, used or unused, COUPLER Device components should be discarded due to compromised sterility. Any packaging or components exposed to human tissue/fluids should be disposed of per hospitals protocols.

INSTRUCTIONS FOR COMPLETION OF IMPLANT CARD:

Enter the following information into the spaces indicated on the Implant Card and then give to patient.

- Name and address of the implanting healthcare institution or provider
- Patient name or patient ID
- Date of implantation

The diagram shows a dashed-line border representing the Implant Card. On the left side, there are several icons and labels: a computer monitor icon with the URL <https://edocs.baxter.com>, a person icon with a plus sign, a person icon with a question mark, a '31' icon, a 'REF' label, a 'LOT' label, and a 'UDI' label. The main area of the card is divided into three horizontal sections by dashed lines, intended for handwritten information.

DISCLAIMER OF WARRANTY:

Synovis Micro Companies Alliance, Inc., a subsidiary of Baxter International Inc. warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. Since SMCA has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SMCA does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SMCA will replace any device which is defective at the time of shipment. No representative of SMCA may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

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Fax : +1 205.941.1522

Site Web : synovismicro.com

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COMPOSANT	REF
COUPLER, 1,0 MM, 1PK	GEM2750/I
COUPLER, 1,0 MM, 6PK	GEM2750
COUPLER, 1,5 MM, 1PK	GEM2751/I
COUPLER, 1,5 MM, 6PK	GEM2751
COUPLER, 2,0 MM, 1PK	GEM2752/I
COUPLER, 2,0 MM, 6PK	GEM2752
COUPLER, 2,5 MM, 1PK	GEM2753/I
COUPLER, 2,5 MM, 6PK	GEM2753
COUPLER, 3,0 MM, 1PK	GEM2754/I
COUPLER, 3,0 MM, 6PK	GEM2754
COUPLER, 3,5 MM, 1PK	GEM2755/I
COUPLER, 3,5 MM, 6PK	GEM2755
COUPLER, 4,0 MM, 1PK	GEM2756/I
COUPLER, 4,0 MM, 6PK	GEM2756

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The diagram shows a dashed-line border representing the implant card. At the top left is a logo and the URL <https://edocs.baxter.com>. Below the URL are three large empty rectangular boxes for text entry. On the left side, there are four small icons: a person with a plus sign, a person with a question mark, a calendar icon with the number 31, and a box labeled REF. Below these are two more boxes labeled LOT and UDI.

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The diagram shows a dashed-line border representing a patient card. At the top left is a small icon of a person with a plus sign. To its right is the URL <https://edocs.baxter.com> and a small icon of an open book. Below these are four rows of input fields, each with a label to its left: 1. A large empty box for patient name and address, with a label 'P.A.' (Patient Address) to its left. 2. A smaller empty box for patient name or ID, with a label 'P?' (Patient Name/ID) to its left. 3. A box for the reference number, with a label 'REF' to its left. 4. A box for the lot number, with a label 'LOT' to its left. At the bottom is a box for the Unique Device Identifier, with a label 'UDI' to its left.

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The diagram shows a dashed-line border representing the implant card. On the left side, there are several icons and labels: a person icon with a plus sign, a person icon with a question mark, a calendar icon with the number 31, a box labeled 'REF', a box labeled 'LOT', and a box labeled 'UDI'. At the top, there is a URL <https://edocs.baxter.com> and a document icon. The main area of the card is divided into three horizontal sections, each with a dashed-line border, intended for text entry.

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The diagram shows a dashed-line border representing the implant card. At the top left is a computer icon with the URL <https://edocs.baxter.com>. At the top right is a document icon. On the left side, there are four icons: a person with a plus sign, a person with a question mark, a box with '31', and a box with 'REF'. Below these are three boxes labeled 'LOT' and 'UDI'. The main area of the card is divided into three horizontal sections by dashed lines, intended for handwritten information.

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The diagram shows a dashed-line border representing the implant card. At the top left, there is a small icon of a person with a plus sign. To its right is the URL <https://edocs.baxter.com>. Below the person icon is a small icon of a person with a question mark. To the right of this icon is a small icon of a document with a question mark. Below these icons are three large empty rectangular boxes for text entry. On the left side of the card, there are four small rectangular boxes labeled 'REF', 'LOT', and 'UDI' (with a small '3' above it). At the top right, there is a small icon of a document with a question mark.

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The diagram shows a dashed-line rectangular label with the following elements:

- Top left: A small icon of a document and the URL <https://edocs.baxter.com>.
- Top right: A small icon of an information book.
- Left side: A vertical column of icons: a person with a plus sign, a question mark, a box with '31', a box with 'REF', a box with 'LOT', and a box with 'UDI'.
- Right side: Three large, empty rectangular boxes stacked vertically, intended for handwritten information.

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The diagram shows an implant card with the following fields and labels:

- Top left: A small icon of a person and a document, with the URL <https://edocs.baxter.com> next to it.
- Top right: A small icon of a document with an 'i' inside.
- Second row: A large empty rectangular box.
- Third row: A small icon of a person with a question mark, followed by a horizontal line.
- Fourth row: A small icon of a document with '31' inside, followed by a horizontal line.
- Fifth row: A small box labeled 'REF'.
- Sixth row: A small box labeled 'LOT'.
- Seventh row: A small box labeled 'UDI'.
- Bottom: A large empty rectangular box.

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The diagram shows a dashed-line border representing an implant card. At the top left is a laptop icon with the URL <https://edocs.baxter.com>. To its right is a document icon with an 'i'. Below these are several fields: a person icon with a plus sign, a person icon with a question mark, a '31' icon, a 'REF' box, a 'LOT' box, and a 'UDI' box. The main body of the card is divided into two large empty rectangular areas for text entry.

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The diagram shows a dashed-line border representing the implant card. At the top left is a small icon of a person with a plus sign. To its right is the URL <https://edocs.baxter.com>. At the top right is a small icon of a document with an 'i'. Below these are three large rectangular boxes with dashed borders, intended for handwritten information. On the left side of the card, there are four small icons: a person with a plus sign, a question mark, a box with '31', and a box with 'REF'. At the bottom left, there are three boxes labeled 'REF', 'LOT', and 'UDI'.

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The diagram shows a dashed-line border containing a template for an implant card. At the top left is a computer icon and the URL <https://edocs.baxter.com>. At the top right is an information icon. Below these are three main sections, each with a dashed border: 1. A large empty box for patient information, with a '+' icon and a person icon to its left. 2. A smaller empty box for patient name or ID, with a person icon and a question mark to its left. 3. A large empty box for implant details, with a '31' icon to its left. To the left of this third box are four vertically stacked labels: REF, LOT, and UDI.

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The diagram shows a dashed-line border representing an implant card. On the left side, there are several icons and labels: a computer icon with the URL <https://edocs.baxter.com>, a person icon with a plus sign, a question mark icon, a calendar icon with the number 31, a box labeled REF, a box labeled LOT, and a box labeled UDI. The main area of the card is divided into three horizontal sections by dashed lines, intended for handwritten information.

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REFERENCES:

1. DeLacure M and Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. Journal of Reconstructive Microsurgery 13:8, 1997.
2. Caution needed when performing MRI scans on patients with aneurysm clips. FDA Medical Bulletin Volume 23, Number 2, June 1993.
3. Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report.

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CH REP

Baxter Healthcare SA
Thurgauerstrasse 130
8152 Glattpark (Opfikon)
Switzerland

 Synovis Micro Companies Alliance, Inc. (a subsidiary of Baxter International Inc.),
439 Industrial Lane, Birmingham, AL 35211-4464 USA
(Tel) +1 205.941.0111 (US Toll free) 1.800.510.3318 (Fax) +1 205.941.1522 www.synovismicro.com

COUPLER and the GEM design are registered trademarks of Baxter International Inc., or its subsidiaries.

MS-0742189/B
2025-12-08

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Baxter Deutschland GmbH