# GEM

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2797

## Symbol Definitions

	Country of Origin		Do Not Use If Package is Damaged
	Refer to Instruction manual/booklet		Temperature - Storage - Operate - Transport
	Content		Humidity - Non-Condensing - Including Condensation
	Separate Collection		Warning
<b>Rx ONLY</b>	CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.		Manufacturer
	Unique Device Identifier		Manufacturer Part Number
	Non-sterile		Authorized Representative in the European Community
	Date of Manufacture with Country of Origin		Authorized Representative in United Kingdom Only
	Catalog number		Authorized Representative in Switzerland Only
	Lot number		
	Medical Device		

## Description

Synovis reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific Instructions provided with each device. Electrosurgical Cords are designed to conduct electrical power from the output of a high frequency generator to a bipolar instrument. Each reusable device is designed to withstand 20 sterilization cycles when properly cared for and sterilized according to the instructions provided below.

## Intended Use

Bipolar Forceps are intended for use in coagulating tissue during surgical procedures and are intended for use with a maximum voltage of 500 Volts (Vpeak).

The function of the cable is to connect the Electrosurgical Forceps directly to the ESU. Bipolar Cables are intended for use with Bipolar Forceps for coagulating tissue during surgical procedures 500 Volts (Vpeak).

## Indications

Indicated to facilitate tissue grasping and control the bleeding of small vessels by bipolar coagulation during surgery.

## Contraindications

None.

## Clinical Benefits

The primary patient benefits are to facilitate a surgical procedure through providing a means for coagulation.

## Patient Population

- Patient should not be user
- Age: Infant to geriatric
- Weight: No restriction

## Intended User

- Education – Trained physician, physician's assistant, nurse, nurse practitioner. No maximum.
- Knowledge: no maximum
  - Understands electrosurgery and electrosurgical techniques.
  - Read and understands supplied User's Guide (Accompanying Document).
  - Understands hygiene.
- Maximum: There is no maximum
- Experience:
  - Minimum:
    - Some training on techniques or training under surveillance/supervision.
    - No special experience needed.
  - Maximum: There is no maximum
- Permissible Impairments:
  - Mild reading / vision impairment or vision correction to 20/20.
  - Impaired by 40% resulting in 60% of normal hearing at 500 Hz to 2.0 kHz.

## Medical Conditions / Disease State

Patients undergoing surgery who will require the use of one or more electrosurgical instruments.

## Handling

- Device use and application to be determined by the physician.
- Remove and discard the tip protector on all forceps prior to use.
- Inspect all devices before each use and send for repair or discard if damaged.

- After the procedure, the cords should be grasped by the plug to disconnect it from the electrosurgical generator. Pulling on the cord to detach it from the generator can cause weaknesses that could lead to a fire hazard.

## Reprocessing Instructions

- Remove excess biological soil from the instruments with a disposable wipe. Place the devices in a container of distilled water or cover with a damp towel.
- Note: Soaking in an enzymatic solution prepared according to the manufacture will facilitate cleaning, especially in instruments with complex features such as lumens, mating surfaces, blind holes, and cannulas
- If instruments cannot be soaked or maintained damp then they should be cleaned as soon as possible (within 60 minutes is recommended) after use to minimize the potential for drying prior to cleaning.
- Used instruments must be transported to the decontamination area for reprocessing in a closed or covered container to prevent unnecessary contamination risk.
- Preparation for cleaning: All cleaning solutions should be prepared at the dilution temperature recommended by the manufacturer. Softened tap water may be used to prepare cleaning solution.
- Note: Fresh cleaning solutions should be prepared when existing solutions become grossly contaminated (turbid).

## Manual Cleaning Steps

- Prepare a proteolytic enzyme solution according to the manufacturers instructions
- Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution to all surfaces. Lumens, blind holes, and cannulas should be flushed with a syringe to remove bubbles and ensure contact of the solution with all

instrument surfaces.

- Soak instruments for a minimum of 10 minutes. While soaking scrub the surfaces using a soft nylon bristle brush until all visible soil has been removed.
- All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.
- Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable, hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
- Prepare and ultrasonic cleaning bath with detergent according to the manufacturer's recommendations. Completely submerge instruments in the cleaning solution and gently shake to them to remove any trapped bubbles. Lumens, blind holes, and, cannulations should be flushed with a syringe to remove bubbles and ensure contact of solution with all instrument surfaces.
- Sonically clean the instruments at the time, temperature, and frequency recommended by the equipment manufacturer and optimal for the detergent used. A minimum of ten (10) minutes is recommended.
- Cords should be cleaned in the ultrasonic cleaner.
- Remove the instruments from the ultrasonic bath and rinse in purified water for a minimum of one(1) minute or until there is no sign of residue detergent or biological soil. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
- Dry instruments with a clean, absorbent, non-shedding, lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas and difficult to access areas.

#### Combination Manual/Automated Cleaning Steps

- Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
- Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulas should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
- Soak instruments for a minimum of 10 minutes. While soaking scrub the surfaces using a soft nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms.
- Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, rough surfaces and areas with moving components or springs.
- Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.
- Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- Place instruments in a suitable validated washer/disinfector. Follow the washer/disinfector manufacturer's instructions for loading the instruments for maximum cleaning exposure, e.g. open all instruments, place concave instruments on their side or upside down, use baskets and trays designed for washers, place heavier instruments on the bottom of the trays and baskets. If the washer/disinfector is equipped with special racks (e.g. for cannulated instruments) use them according to the manufacturer's instructions.
- Process instruments using a standard washer/disinfector instrument cycle

according to the manufacturer's instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-Wash - Cold Softened Tap Water - 2 minutes
2	Enzyme Spray and Soak - Hot Softened Tap Water - 1 minute
3	Rinse - Cold Soften Tap Water
4	Detergent Wash - Hot Tap Water (64-66°C/146-150°F) - 2 minutes
5	Rinse - Hot Purified Water (64-66°C/146-150°F) - 1 minutes
6	Hot Air Dry(116°C/146-240°F) - 7 to 30 minutes

#### The washer/disinfector manufacturer's instructions should be followed.

- A washer/disinfector with demonstrate efficacy (e.g. FDA approval, validated to ISO 15883) should be used.
- Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfector.
- Many manufacturers pre-program their washer/disinfector with standard cycles that may include a thermal low-level disinfection cycle after the detergent wash. The thermal disinfection cycle should be performed to achieve a minimum value A0 = 600 (e.g. 90°C/194°F for 1 minute according to ISO 15883-1) and is compatible with applicable instruments.
- Dry instruments with a clean, absorbent non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

#### Sterility

- Single Devices may be packaged in an approved (e.g. FDA cleared or ISO 11607 compliant) medical grade sterilization pouch or wrap. Care should be used when

packaging so that the pouch or wrap is not torn. Devices should be wrapped using a double wrap or equivalent method (ref: AAMI ST79, AORN Guidelines).

- Reusable wraps are not recommended.
- Instruments may be packed in an approved (e.g. FDA cleared or ISO 11607 compliant) general use perforated tray or case along with other devices under the following conditions.
  - Arrange all devices to allow access of steam to all surfaces. Open hinged devices and ensure devices are disassembled if it is recommended.
  - The case or tray must be wrapped in an approved (e.g. FDA cleared or ISO 11607 compliant) medical grade sterilization wrap by following the double wrap method or equivalent (ref: AAMI ST79, AORN Guidelines).
- Follow the case/tray manufacturer's recommendations for loading and weight. Total weight of wrapped case or tray should not exceed 11.4kg/25lbs.
- Instruments may be packaged in an approved (e.g. FDA cleared or ISO 11607 compliant) rigid container system (i.e. those with filters or valves) along with other devices under the following conditions:
  - The container manufacturer's recommendations should be followed regarding preparation, maintenance and use of the container.
  - Arrange all devices to allow access of steam to all surfaces. Open hinged devices and ensure devices are disassembled if it is recommended.
- Follow the container manufacturer's recommendations for loading and weight. Total weight of a filled container system should not exceed 11.4kg/25lbs.
- Moist heat/steam sterilization is the recommended method for Synovis instruments.
- Use of an approved chemical integrator (class 5) or chemical emulator (class 6) within each sterilization is recommended.

- Always consult and follow the sterilizer manufacturer's instructions for load configuration and equipment operation. Sterilization equipment should have demonstrated efficacy (e.g. FDA clearance, EN 13060 or EN 285 compliance). Additionally, the manufacturer's recommendations for installation, validation, and maintenance should be followed.
- Validated exposure times and temperatures are listed in the following table:

Cycle Type	Minimum Temperature	Minimum Exposure Time
United States Recommended Parameters		
Pre-Vacuum/ Vacuum Pulse	132°C / 270°F	4 Minutes
Cycle Type	Minimum Temperature	Minimum Exposure Time
European Recommended Parameters		
Pre-Vacuum/ Vacuum Pulse	134°C / 273°F	3 Minutes

### Drying and Cooling

- The recommended drying time for single wrapped instruments is 20 minutes unless otherwise noted in device specific instructions.
- Drying times for instruments processed in containers and wrapped trays can vary depending upon the type of packaging, type of instruments, type of sterilizer, and total load. A minimum dry time of 30 minutes is recommended, but to avoid wet packs, extended dry times greater than 30 minutes may be needed for larger loads under certain conditions or if otherwise recommended in accompanying documentation. For large loads verification of dry times by the health care provider is recommended.
- A 30 minute minimum cooling time is recommended after drying but longer times may be necessary because of load configuration, ambient temperature and humidity, device design and packaging used.

### Storage

Reusable bipolar forceps and cords must be stored individually in a protective container in a clean, dry, moisture free area with suitable protection from extreme temperature and humidity. Protect cords from mechanical damage and kinking and handle with care.

### Notices

- Electrosurgical equipment and accessories are intended to be used by health professional educated in their use.
- After Use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.
- Any serious incident that has occurred in relation to the device should be reported to Synovis, the governing body, and the Competent authority of the member state in which the incident occurred.
- Please refer to the manufacturer of the generator for warnings, precautions, contraindications, undesirable side-effect, measures to be taken, and limitations of use for the electrosurgical system and accessories
- Reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific Instructions provided with each device.

### Cautions

- Active electrosurgical devices must be placed away from the patient when not in use. A holster is recommended.
- Use a wet sponge to remove eschar build up from the tips of the forceps.
- Consult the RF generator Owner's Manual. Use only in the active RF generator bipolar output connections.

- Do not perform electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids, or objects, or in the presence of oxidizing agents, as a fire could result.
- Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.
- Do not bend, cut, or modify any electrosurgical device before or during surgery.
- Keep the voltage/power settings at the lowest possible level to achieve the desired effect.
- Establish a protocol for device examination and maintenance.
- Electromagnetic compatibility (EMC) – Aspen Surgical devices are considered medical electrical equipment. Medical electrical equipment may require precautions regarding EMC and may need to be installed or put into service according to the EMC information provided in the operator's manual of the generator.
- Do not use electrosurgical equipment unless properly trained to use it in the specific procedure being undertaken. Use by physicians without such training may result in serious, unintended patient injury.

### Warning

- Improper use of any element in electrosurgical systems may cause serious burn. The use of the bipolar forceps with safely functioning and compatible cords, electrosurgical generators and accessories is the user's responsibility. Read and understand all precautions and instructions before using any electrosurgical device.
- Avoid HF output mode/power settings where maximum output voltage may exceed the stated maximum voltage.

- Do not make any modification to these devices. Modification may cause burn or injury.
- Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols which may include toxic gases and vapors, live and dead cellular material, and viruses), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures. Contact Aspen Surgical at [customerservice@aspensurgical.com](mailto:customerservice@aspensurgical.com) for additional information or to inquire about our smoke evacuation solutions.

### Warranty

All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration to any product that is performed by anyone not authorized by Synovis will result in the immediate loss of warranty.

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Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

**CH REP**

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

**UK REP**

Emergo Consulting (UK) Limited  
c/o Cr360 – UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ,  
United Kingdom

Manufactured for: Synovis Micro Companies Alliance, Inc. (A subsidiary of Baxter International Inc.)  
439 Industrial Lane, Birmingham, AL 35211-4464 USA  
205.941.0111 (USA) • 800.510.3318 • [www.synovismicro.com](http://www.synovismicro.com)

 Aspen Surgical Products Inc., 6945 Southbelt Dr SE, Caledonia, MI 49316 USA  
Phone +1 616-698-7100 | Toll-Free 888-364-7004 | [aspensurgical.com](http://aspensurgical.com)

0001331 REV. 18  
2025-08-11